

Gryllus Techniques and Protocols**in situ Hybridization Probe Synthesis**

To generate probes for use visualizing mRNA in Gryllus.

General Notes:

- 1) Amplify your target gene.
- 2) Assemble the following reaction in a 1.5mL eppendorf tube, using your target gene as template and either T7 or Sp6-T7 primers for sense and antisense probes. Seal tube lids with parafilm.

Reagent	1x (μl)
10x RNA polymerase buffer	2
Dig labeling mix	2
RNase inhibitor	1
RNA Polymerase	2
500ng template	x
MQ H ₂ O	To 20μl

- 3) Incubate 3-4h at 37°C.
- 4) Add 30μl TE, 5μl 3M NaOAc, and 2.5volumes (125μl) of cold 100%EtOH.
- 5) Incubate on ice at least 20min.
- 6) Spin 15min at max speed at 4°C in a microfuge.
- 7) Aspirate and discard supernatant.
- 8) Add 500μl; cold 75% EtOH in H₂O.
- 9) Invert several times to mix.
- 10) Spin 5min at max speed at 4°C in a microfuge.
- 11) Aspirate and discard supernatant.
- 12) Air dry 2-3 (and certainly under 5) minutes.
- 13) Dissolve in 20μl TE. Do not pipette or vortex.
- 14) Dilute 1:10 in TE to check reaction by measuring the concentration on a spectrophotometer and running in gel electrophoresis.
*Wash box and comb, and change running buffer before making gel to decrease RNases.
- 15) Adjust to 100ng/μl in hybridization solution. Store at -20°C.