

January 2014

RNAi Handbook

From dsRNA synthesis
to RNAi Treatment



NSF/EDEN Grant number IOS # 0955517

Required Equipment

Sample	
Attached gemmules in a 12 wells plate	
Solutions	
Strekal's medium	1mL/well/day
Strekal's medium (for wash)	2mL/well
Primers stock solutions (100µM)	2µL/reaction
Sponge cDNA	1µL/gene
3M Sodium Acetate (pH 5.2)	2.2µL/reaction
Isopropanol	22µL/reaction
RNAlater® solution	200µL/treatment
Kits	
GoTaq® Green Master Mix (or any similar Kit)	
Cloning Kit	
QIAquick® Gel Extraction Kit (or any similar Kit)	
QIAquick® PCR purification kit (or any similar Kit)	
QIAprep® Spin Miniprep Kit (or any similar Kit)	
Promega T7 RiboMAX™ Express RNAi System	
Qiagen RNeasy® Mini Kit (or any similar Kit)	
Invitrogen Super Script® III Reverse Transcriptase	
Tools	
Thermocycler	
Incubator	

Primers design

Several primers will be required to synthesize the dsRNA and check the knockdown by qPCR :

- One pair to target the 5' end of the mRNA (product size ~200pb)
- Primers as above with a partial T7 sequence (5'CGACTCACTATAGGG...3')
- Full T7 primer (5'ATAGAATTCTCTAGAAGCTTAATACGACTCACTATAGGG3')
- One primer pair targeting the 3' end of the mRNA (product size ~100pb)

dsRNA synthesis

1. Perform a regular PCR using primers targeting 5' end of your gene of interest (PCR 0).

For example using GoTaq® Green Master Mix

Mix	
GoTaq	25µL
Fw-Primer	2µL
Rv-Primer	2µL
cDNA (~200ng)	1µL
H2O	20µL

Cycles	
95°C	5'
95°C	15s
55°C	30s
72°C	30s
72°C	5'

) x30

2. Check 2µL of PCR product on a 1% agarose gel.
3. Purify DNA using QIAquick® PCR purification kit.
Or running all the PCR product and purify it with QIAquick Gel Extraction Kit depending on gel profile.
4. Proceed with a plasmid ligation (e.g., TOPO cloning) and a bacterial transformation.
5. After checking colonies by PCR, carry out an overnight cell culture and extract plasmid with QIAprep® Spin Miniprep Kit.
Sequence it in order to confirm the insert.
6. Perform a PCR as in step 1 with primers containing the partial T7 extension and the miniprep solution (or the PCR 0 product) as template (PCR 1).
7. As previously, run 2µL of PCR product and purify DNA.
8. Perform a PCR as in step 1 with Full T7 primers and DNA from step 7 as template (PCR 2).
9. One more time, check 2µL of PCR product and purify DNA.

10. Quantify DNA concentration of PCR 2.

Need about 1 μ g for dsRNA synthesis.

11. Prepare the T7 RiboMAX™ mix as following:

Mix	
RiboMAX™	10 μ L
PCR 2 product	1 μ g
Nuclease-Free Water	x μ L
Enzyme Mix, T7 Express	2 μ L
Final Volume	20 μ L

12. Incubate at 37°C for 30 minutes.

13. Remove DNA template by adding 1 μ L freshly diluted Rnase Solution and 1 μ L RQ1 Rnase-Free Dnase. Then, incubate 30 minutes at 37°C.

According to the protocols supplied by Promega.

14. Purify dsRNA with 2.2 μ L of 3M Sodium Acetate (pH 5.2) and 22 μ L of isopropanol.

15. Incubate overnight at -20°C

16. Spin at top speed for 10 minutes.


17. Remove supernatant and wash pellet with cold 70% ethanol.

18. Remove ethanol, air dry the pellet for 15 minutes and resuspend dsRNA in 50 μ L of Nuclease-Free Water.

19. Quantify dsRNA solution and run 2 μ L of both dsRNA and PCR 2 products.

A shift should appear between the PCR 2 product and dsRNA migrations (see picture 1).

20. Aliquot and store at -20°C.



RNAi Treatment

1. Check that gemmules are at the desire stage of development and take pictures.
2. Dilute dsRNA in 1mL of new Strekal's medium (10µg/mL final concentration).
3. Well by well, quickly but gently replace de medium with the medium containing the dsRNA.
4. Change all the solutions every 24h.
Also in control wells.
5. Take pictures as often as required.
6. Remove dsRNA after 48h and wash sponges two times with 1 mL of Strekal's medium.
7. Using a 1mL pipette, collect sponges in a 1.5mL microcentrifuge tubes combining like treatments in the same tube.

- 8. Spin 3 minutes at top speed, remove liquid and add 200 μ L of RNeasy[®] and store at -20 $^{\circ}$ C until used.**
- 9. Extract mRNA using Qiagen RNeasy[®] Mini Kit with a step of DNA degradation and quantify.**
- 10. Synthesize corresponding cDNA from 125ng of mRNA as template using Invitrogen Super Script[®] III Reverse Transcriptase.**
- 11. Check gene expression shutdown by qPCR using primers targeting the 3' end of the mRNA.**

April Hill, *Department of Biology, University of Richmond, Richmond, VA, USA*
(ahill2@richmond.edu)

Quentin Schenkelaars, *Institut Méditerranéen de Biodiversité et d'Ecologie marine et continentale (IMBE), Aix Marseille-Université, Marseille, France* (quentin.schenkelaars@univ-amu.fr)

Fernando Gomez, *Department of Biology, University of Richmond, Richmond, VA, USA*
(fgomez@richmond.edu)

Chelsea Hall, *Department of Biology, University of Richmond, Richmond, VA, USA*
(chelsea.hall@richmond.edu)

Laura Fierro, *Institut Méditerranéen de Biodiversité et d'Ecologie marine et continentale (IMBE), Aix Marseille-Université, Marseille, France* (laura.fierro@univ-amu.fr)

Emmanuelle Renard, *Institut Méditerranéen de Biodiversité et d'Ecologie marine et continentale (IMBE), Aix Marseille-Université, Marseille, France* (emmanuelle.renard@univ-amu.fr)

Carole Borchiellini, *Institut Méditerranéen de Biodiversité et d'Ecologie marine et continentale (IMBE), Aix Marseille-Université, Marseille, France* (carole.borchiellini@univ-amu.fr)



NSF/EDEN Grant number IOS # 0955517