

## Protocols:

Trainee: Anna Vickrey

Host lab: Shapiro Lab, University of Utah

Dates of Internship: 06/26/2012- 09/1/2012

Title: Gibson Cloning (from Dr. Eric Domyan)

Rationale and background: To ligate insert and vector fragments for high efficiency plasmid transformations.

Protocol:

I: Quick Protocol

1. Thaw mix on ice.
2. Split into two tubes (7.5ul each)
3. Add 2.5ul total DNA (target 1:1:1... molar ratios but this is not critical)
4. Incubate at 50C for 1hr
5. Transform 2 ul in 50-100 ul of cells.

\*If you want to chew up methylated background DNA add .5 ul DpnI and Incubate at 37 for 15 -30 min

II: Gibson reaction from set up:

1. Prepare 5X ISO buffer. Six ml of this buffer can be prepared by combining the following:

- 3 ml of 1 M Tris-HCl pH 7.5
- 150  $\mu$ l of 2 M MgCl<sub>2</sub>
- 60  $\mu$ l of 100 mM dGTP
- 60  $\mu$ l of 100 mM dATP
- 60  $\mu$ l of 100 mM dTTP
- 60  $\mu$ l of 100 mM dCTP
- 300  $\mu$ l of 1 M DTT
- 1.5 g PEG-8000
- 300  $\mu$ l of 100 mM NAD
- Add water to 6 ml

2. Aliquot 100  $\mu$ l and store at -20 °C

3. Prepare an assembly master mixture. This can be prepared by combining the following:

- 320  $\mu$ l 5X ISO buffer
- 0.64  $\mu$ l of 10 U/  $\mu$ l T5 exonuclease
- 20  $\mu$ l of 2 U/ $\mu$ l Phusion polymerase
- 160  $\mu$ l of 40 U/ $\mu$ l Taq ligase
- Add water to 1.2 ml

4. Aliquot 15  $\mu$ l and store at -20 °C. This assembly mixture can be stored at -20 °C for at least one year.

The enzymes remain active following at least 10 freeze-thaw cycles.

This is ideal for the assembly of DNA molecules with 20-150 bp overlaps. For DNA molecules overlapping by larger than 150 bp, prepare the assembly mixture by using 3.2  $\mu$ l of 10 U/  $\mu$ l T5 exo.

5. Thaw a 15  $\mu$ l assembly mixture aliquot and keep on ice until ready to be used.

6. Add 5  $\mu$ l of DNA to be assembled to the master mixture. The DNA should be in equimolar amounts. Use 10-100 ng of each ~6 kb DNA fragment. For larger DNA segments, increasingly proportionate amounts of DNA should be added (e.g. 250 ng of each 150 kb DNA segment).

7. Incubate at 50 °C for 15 to 60 min (60 min is optimal).

8. If cloning is desired, transform 2  $\mu$ l of the reaction mixture into 50-100  $\mu$ l competent cells (DH5 $\alpha$ ) by heat shock