

Protocols:

Trainee: Anna Vickrey

Host lab: Shapiro Lab, University of Utah

Dates of Internship: 06/26/2012- 09/1/2012

Title: Gibson Cloning (from Dr. Eric Domyan)

Rationale and background: To ligate insert and vector fragments for high efficiency plasmid transformations.

Protocol:

I: Quick Protocol

1. Thaw mix on ice.
2. Split into two tubes (7.5ul each)
3. Add 2.5ul total DNA (target 1:1:1... molar ratios but this is not critical)
4. Incubate at 50C for 1hr
5. Transform 2 ul in 50-100 ul of cells.

*If you want to chew up methylated background DNA add .5 ul DpnI and Incubate at 37 for 15 -30 min

II: Gibson reaction from set up:

1. Prepare 5X ISO buffer. Six ml of this buffer can be prepared by combining the following:

- 3 ml of 1 M Tris-HCl pH 7.5
- 150 μ l of 2 M MgCl₂
- 60 μ l of 100 mM dGTP
- 60 μ l of 100 mM dATP
- 60 μ l of 100 mM dTTP
- 60 μ l of 100 mM dCTP
- 300 μ l of 1 M DTT
- 1.5 g PEG-8000
- 300 μ l of 100 mM NAD
- Add water to 6 ml

2. Aliquot 100 μ l and store at -20 °C

3. Prepare an assembly master mixture. This can be prepared by combining the following:

- 320 μ l 5X ISO buffer
- 0.64 μ l of 10 U/ μ l T5 exonuclease
- 20 μ l of 2 U/ μ l Phusion polymerase
- 160 μ l of 40 U/ μ l Taq ligase
- Add water to 1.2 ml

4. Aliquot 15 μ l and store at -20 °C. This assembly mixture can be stored at -20 °C for at least one year.

The enzymes remain active following at least 10 freeze-thaw cycles.

This is ideal for the assembly of DNA molecules with 20-150 bp overlaps. For DNA molecules overlapping by larger than 150 bp, prepare the assembly mixture by using 3.2 μ l of 10 U/ μ l T5 exo.

5. Thaw a 15 μ l assembly mixture aliquot and keep on ice until ready to be used.

6. Add 5 μ l of DNA to be assembled to the master mixture. The DNA should be in equimolar amounts. Use 10-100 ng of each ~6 kb DNA fragment. For larger DNA segments, increasingly proportionate amounts of DNA should be added (e.g. 250 ng of each 150 kb DNA segment).

7. Incubate at 50 °C for 15 to 60 min (60 min is optimal).

8. If cloning is desired, transform 2 μ l of the reaction mixture into 50-100 μ l competent cells (DH5 α) by heat shock

Title: RNA in situ hybridization (from Dr. Eric Domyan)

Rationale and background: To visualize gene expression by binding a tagged probe to target mRNA

Protocol:

DAY 0 – Dissect tissue, fix in 4% paraformaldehyde (4% PFA) at 4°C overnight. Dehydrate through methanol gradient into 100% MeOH at 4°C, store at -20°C.

DAY 1 – All incubation steps to be performed with gentle agitation on an orbital shaker, unless otherwise indicated.

1. Prepare supplies and preheat solutions.
 - A. Transfer prehyb stock solution into sealed conical tube and preheat to 60.5°C.
2. Rehydrate fixed tissue (In 6-well baskets or nunc tubes. Use 5mL/wash for baskets and 2mL/wash for nunc tubes.)
 - A. Wash tissue 1 x 15 min each in 75%MeOH/25% PBSTw, 50%MeOH/50%PBSTw, 25% MeOH/75% PBSTw.
 - B. Wash tissue 2 x 15 min in PBSTw.
3. Incubate tissues for 30 min at 25°C in 2mL/well of 6% H₂O₂.
4. Wash tissues 4 x 5 min at 25°C with 2mL/well PBSTw.
5. Incubate tissues 12 min at 25°C in 2mL/well proteinase K solution.
6. Wash tissues 1 x 5 min at 25°C with 2mL/well PBSTw.
7. Incubate tissues 20 min at 25°C in 2mL/well post-fix solution.
8. Wash tissues 2 x 5 min at 25°C with 2mL/well PBSTw.
9. Transfer tissues into separate 2mL Nunc tubes, 1 tube per probe.
10. Prehybridization step: Incubate for at least 1 hr at 60.5°C in 1.5 to 2mL prehyb buffer.
11. Hybridization step: Add 0.65mg probe/well and incubate sections overnight in hybridization chamber at 60.5°C in probe+prehyb buffer.
12. Prepare and preheat solutions for Day 2: Add SDS to solution 1 (to a final concentration of 1%) and preheat overnight at 60.5°C.

DAY 2 - All incubation steps to be performed with gentle agitation on an orbital shaker, unless otherwise indicated.

1. Transfer tissue into baskets of a 6 well plate. Remove hyb buffer from each well and wash tissues 3 X 30 min at 60.5°C with 5mL/well Solution 1.
2. Prepare and preheat solutions.
 - A. Prepare 50/50% mix of Solution1/Solution 2 and preheat to 60.5°C.
 - B. Add RNase to Solution 2 and preheat to 37°C.
 - C. Preheat one wash volume of Solution 3 to 25°C and two wash volumes to 60.5°C.
3. Wash tissues 1 X 10 min at 60.5°C with 5mL/well of 50/50% mixture of Solution 1/ Solution 2.
4. Wash tissues 4 X 10 min at 25°C with 5mL/well Solution 2.
5. Incubate tissues 15 min at 37°C in 5mL/well RNase solution.
6. Wash tissues 1 X 10 min at 25°C with 5mL/well Solution 2.
7. Wash tissues 1 X 10 min at 25°C with 5mL/well Solution 3.
8. Wash tissues 2 X 1 hr at 60.5°C with 5mL/well Solution 3.
9. Preheat Tissue Blocking (TB) buffer to 25°C.
10. Wash tissues 3 X 10 min at 25°C with 2ml/well TBSTw.
11. Transfer tissues into new 2mL Nunc tubes.
11. Blocking step
 - A. Incubate tissues at least 2 hr at 25°C in 2mL/well TB buffer.
12. Antibody step
 - A. Dilute anti-DIG AP antibody 1:2000 in TB buffer.

- **B.** Add 2 mL of TB containing the antibody to each Nunc tube.
- **C.** Incubate tissues overnight at 4°C in 2mL TB buffer + antibody.

Day 3 - All incubation steps to be performed with gentle agitation on an orbital shaker, unless otherwise indicated.

1. Remove antibody solution and save (can be reused).
 - **A.** Add sodium azide to a final concentration of 0.2mM to prevent microbial growth.
 - **B.** Store antibody solution at 4°C and reuse up to two additional times.
2. Transfer tissues into baskets of 6 well plate.
3. Wash tissues 8 X 30 min at 25°C with 5mL/well TBSTw containing 2 mM levamisole.
4. Remove tissues from baskets, use forceps to separate sections & remove visible debris, transfer tissues into clean Nunc tubes.
5. Wash tissues 1 X 10 min at 25°C with 1.5mL NTMT containing 2 mM levamisole.
6. Detection step
 - **A.** Protect tissues from light and incubate at 25°C in 1mL/tube of a 60/40% mixture of NTMT containing 2mM levamisole/BM Purple.
 - **B.** Color development time ranges from several hours to several days. Change NTMT/BM Purple solution as needed (substrate will precipitate over time).
 - Once color is fully developed, wash tissues 2 X 5 min at 25°C with 1mL/tube NTMT containing 2mM levamisole.
7. Bleaching step
 - **A.** Post-fix tissues overnight at 4°C in 1mL/tube 4% PFA.
 - **B.** Remove PFA and incubate tissues for 30 min at 25°C in 1mL/tube PBSTw containing 3% H₂O₂.
 - **C.** Wash tissues 1 X 10 min at 25°C in 1mL/tube PBSTw.
 - **D.** Store tissues at 4°C in 1mL/tube 4% PFA.

Solutions

Day 1

PBSTw: 1X PBS containing 0.1% Tween-20 and 0.2mM sodium azide. Pass through 0.2µm filter to remove insoluble material/contaminants.

6% H₂O₂: 1mL 30% H₂O₂ per 4mL PBS

Proteinase K: 0.25µL 20mg/mL proteinase K per 1mL PBSTw

Post-fix: 8µL 25% glutaraldehyde per 1mL 4% PFA

Prehybridization solution

stock solution	final conc	for 1000mL
100% formamide	50%	500mL
20X SSC	5X	250mL
Blocking reagent	1%	10g
10mg/mL yeast tRNA	10µg/mL	1mL
10mg/mL heparin	10µg/mL	1mL
dH ₂ O to vol.	--	to 1000mL

Aliquot 50mL volumes into conical tubes & store at -20°C.

Solution 1

stock solution	final conc	for 500mL
----------------	------------	-----------

100% formamide	50%	250mL
20X SSC	5X	125mL
dH ₂ O	--	75mL
10% SDS	1%	50mL

Mix formamide, SSC & dH₂O as shown above & store at –20°C.

DO NOT add SDS to freezer stock solution because SDS will precipitate in the cold.

Prior to using Solution 1 for ISH, add 1mL of 10% SDS per 9mL Solution 1 stock.

Day 2

Solution 2

stock solution	final conc	for 500mL
1M Tris-HCl, pH 7.5	10mM	5mL
5M NaCl	0.5M	50mL
100% Tween-20	0.1%	0.5mL
0.2M sodium azide	0.2mM	0.5mL
dH ₂ O	--	444mL

Pass through 0.2µm filter to remove insolubles/contaminants, store at 25°C.

RNase: 5µL RNase (50µg/mL) per 1mL Solution 2

Solution 3

stock solution	final conc	for 500mL
20X SSC	2X	50mL
100% formamide	50%	250mL
dH ₂ O	--	200mL

Store at –20°C

Tissue blocking (TB) buffer

stock solution	final conc	for 500mL
10X TBS	1X	50mL
100% sheep serum	10%	50mL
10% blocking reagent	1%	50mL
BSA	1%	0.5g
dH ₂ O to vol.	--	to 500mL
100% Tween-20	0.1%	0.5mL

Mix TBS, serum, blocking reagent, BSA & dH₂O as shown above.

Filter through #2 Whatman filters.

Add Tween-20.

Aliquot 6mL volumes into conical tubes & store at –20°C.

Sheep serum (must be heat-inactivated before use).

To heat inactivate:

- thaw new bottle of serum
- incubate at 70°C for 30 min
- aliquot & store at -20°C

10% Blocking reagent

stock solution	final conc	for 100mL
maleic acid	100mM	1.2g
5M NaCl	150mM	3mL
dH ₂ O to vol.	--	to 100mL
Blocking reagent	10%	10g

Mix maleic acid, NaCl & dH₂O according to above, pH to 7.5 (note: strong buffer so difficult to pH, try using solid NaOH pellets to raise pH initially).

Add blocking reagent, microwave briefly to aid solubility (avoid boiling over, solution will be cloudy & viscous so watch carefully to ensure blocking reagent is completely dissolved in solution).

Aliquot 10 mL volumes into conical tubes & store at -20°C.

Day 3

2M Levamisole: Dissolve 4.82g levamisole in ~7mL double-distilled H₂O (total volume should equal 10mL), aliquot 200µL volumes & store stocks at -20°C.

TBSTw + levamisole: 1X TBSTw containing 2mM levamisole.

NTMT + levamisole (inhibits endogenous alkaline phosphatases):

stock solution	final conc	for 500mL
1M Tris-HCl, pH 9.5	100mM	50 mL
5M NaCl	100mM	10mL
1M MgCl ₂	50mM	25mL
dH ₂ O	--	415mL
100% Tween-20	0.1%	0.5 mL
2M Levamisole	2mM	0.5mL

Mix Tris, NaCl, MgCl₂ & dH₂O as shown above.

Pass through 0.2µm filter to remove insoluble material/contaminants. Store at 25°C.

DO NOT add Tween or levamisole to stock solution.

Prior to using NTMT for ISH, add 1µL of 100% Tween-20 and 1µL of 2M levamisole per 1mL NTMT stock.

3% H₂O₂: 1mL 30% H₂O₂ per 9mL PBSTw.

Abbreviations

BSA = bovine serum albumin

DIG = digoxigenin

H₂O₂ = hydrogen peroxide

PFA = paraformaldehyde

PBS = phosphate-buffered saline

PBSTw = 1X PBS + 0.1% Tween-20

SDS = sodium dodecyl sulfate (aka lauryl sulfate)

SSC = saline sodium citrate

TBS = Tris-buffered saline

TBSTw = 1X TBS + 0.1% Tween-20

Reagents and Supplies

Anti-DIG antibody, Fab fragments, cat # 11214667001, Roche

Blocking reagent, cat # 11096176001, Roche

BM Purple AP substrate, precipitating, cat # 11442074001, Roche

BSA, cat # BP1600-100, Fisher

Formamide, cat # F5786-1L, Sigma

Glutaraldehyde, 25% solution in H₂O, cat # G6257-100ML, Sigma

Heparin, sodium salt, cat # H3393, Sigma

Hydrogen peroxide, 30% solution in H₂O, cat # BP2633-500, Fisher

Levamisole, cat # L9756, Sigma

Magnesium chloride, cat # M33-500, Fisher

Maleic acid, cat # M0375-500G, Sigma

Paraformaldehyde, cat # 101176-014, VWR

PBS, w/out Ca & Mg, MP Biomedicals powdered media, cat # ICN1760420, Fisher

Proteinase K solution, 20mg/mL, biotechnology grade, cat # E195-5ML, Amresco

RNase, cat # R6513, Sigma

SDS, cat # S529-500, Fisher

Sheep serum, cat # S2263-500mL, Sigma

Sodium azide, granular, cat # S2271-100, Fisher

Sodium chloride, cat # BP358-212, Fisher

SSC, 20X solution, cat # S24022-4000.0, Research Products International

Tris-HCl, cat # BP153-1, Fisher

Tween-20, cat # BP337-100, Fisher

Yeast tRNA, cat # 109495, Roche

Polyester mesh, 33 micron, 12" x 24", cat # CMY-0033-D, Small Parts Inc.

Title: Bacterial Toxicity Assay (from Dr. Eric Domyan)

Rationale and background: Test the toxicity (to *E. coli*) of IPTG-induced expression of different plasmid inserts.

Protocol:

1. For each construct you want to test, add 15 mL LB + 15uL 100mg/mL carbenicillin to a 15mL culture tube.
2. Pick a single colony from a plate with a micropipette tip; swirl tip in your tube of LB+ carb.
3. Incubate tube of LB at 37°C with shaking at 300rpm for 1 hr.
4. Split the culture into 6 x 2 mL aliquots in 15 mL culture tubes.
5. To 3 of the tubes, add 45uL 100mg/mL IPTG.
6. Incubate tubes O/N with shaking at 37°C.
7. The next day, determine the OD600 with a spectrophotometer.

LB recipe: pH=7.0 (1L of ddH₂O)

10g NaCl

10g Bacto Tryptone

5g Yeast Extract

1. Fill container with 750L of ddH₂O

2. Weigh out ingredients
3. Use stir bar to dissolve everything into H₂O.
4. pH to 7.0 with 10M NaOH (approx. 1 drop)
5. Bring volume up to 1L
6. Pour solution into 2-500ml bottles
7. Label solution and date
8. Put autoclave tape on and autoclave for 40 min.

Title: Assay for soluble protein expression (from Dr. Martin Horvath)

Rationale and background: Grow cells containing a construct to test for the expression of soluble GST-tagged protein in *E. coli*. Expect some optimization (induction times and temperatures, buffers, cell growth at induction, etc.). Here, we tested protein expression in 60mL cultures under two induction conditions: 20 hours at 20°C and 3 hours at 37°C. You will need to grow separate 60mL cultures for each of your induction conditions. Use BL21 pLysS strain.

Protocol:

1. Start a 2mL liquid culture from a colony containing the GST-construct that you plan to induce expression in. Grow overnight (O/N) shaking at 180 rpm at 37°C in 2mL LB, 3.7uL 100mg/mL carbenicillin (carb), and 2uL 68mg/mL chloramphenicol (chlor).
2. The next morning, wash the cells and use them to inoculate a 60mL 2xYT culture in a **baffled** flask. Before inoculating, set aside an aliquot (~1mL) of 2xYT to use as a spectrophotometer blank.
 - Spin the 2mL O/N culture at 4750rpm at 4°C for 5 min.
 - Pour off the media, and borrow 2mL 2xYT to re-suspend cells.
 - Spin the 2mL culture at 4750rpm at 4°C for 5 min. Pour off media.
 - Re-suspend again in 2mL borrowed 2xYT, transfer to inoculate ~60mL 2xYT culture with the washed cells.
3. Grow ~60mL culture at 37°C with shaking.
4. Monitor the OD at 600nm starting at least 1 hr. after inoculation.
5. When OD=1 (0.85-1.2 will do), at the same time:
 - Induce expression by adding 300uL IPTG.
 - Move culture to desired expression conditions (e.g. 20hr. at 20°C).
6. After induction, spin down 45mL culture by centrifugation for 10 min at 4750 rpm, 4°C.
7. Pour off the supernatant and re-suspend in 2mL Cell Wash with 7.5uL 0.5M PMSF.
8. Divide into three aliquots:
 - 60uL (After freezing, lyse and run on an SDS-PAGE gel to visualize crude extract)
 - 2X 0.97mL (After freezing, sonicate and purify GST-tagged protein by Batch method, then run on an SDS-PAGE gel to visualize soluble GST-tagged protein expression).
9. Spin down aliquots by centrifugation at room temperature at 6000 rpm.
10. Remove supernatant and freeze at -20°C O/N.

Solutions:

1. LB: see "Toxicity Assay"
2. 2xYT media:
 1. Measure ~900ml of distilled H₂O.
 2. Add 16g Bacto Tryptone.
 3. Add 10g Bacto Yeast Extract.
 4. Add 5g NaCl.
 5. Adjust pH to 7.0 with 5N NaOH.

6. Adjust to 1L with distilled H₂O.

Sterilize by autoclaving.

On the same day that you start a 60mL culture, add 1.8 mL 1 M potassium phosphate pH 7.8, 300uL 1M glucose, 120uL 100mg/mL carb, 60uL 68mg/mL chlor to 60mL 2xYT media.

3. Cell Wash (10 mM Tris pH 8.2, 100 mM NaCl, 1 mM EDTA) makes 500 mL

Combine at room temperature; bring to 500mL total volume with di-H₂O:

- Tris base 0.61 g
- 5 N hydrochloric acid (HCl) 0.5 mL
- 0.5 M EDTA 1.0 mL

Do not add PMSF to stock Cell Wash.