

Anna Bigney

Lane Lab, University of Rhode Island

June-August 2014

Red and Brown Algal DNA Extraction

Purpose: To extract algal genomic DNA

Procedure:

1. Grind .5 μ L of sample in liquid nitrogen
2. Add 500 μ L of Red Algal Buffer for Red Algae or 500 μ L of Brown Algal Buffer for Brown Algae
3. Add 10 μ L of 10mg/mL ProK
4. Add 50 μ L 10% Tween
5. Incubate on rotator for an hour
6. Freeze at -20 °C for twenty minutes
7. Centrifuge at 13.2 KxG for ten minutes
8. Warm Elution Buffer (Buffer PE) to 65 °C
9. Pipette supernatant into filtration column
10. Centrifuge at 11 KxG for two minutes
11. Add 450 μ L of Buffer PC to flow-through and pipette up and down five times to mix
12. Pipette 700 μ L supernatant/Buffer PC solution into binding column
13. Centrifuge at 11 KxG for one minute and discard flow-through
14. Repeat previous two steps for remainder of solution
15. Apply 700 μ L Buffer PW2 Wash solution
16. Centrifuge at 11 KxG for one minute and discard flow-through
17. Apply additional 200 μ L of Buffer PW2 Wash solution
18. Centrifuge at 11 KxG for two minutes
19. Place column in fresh 2mL tube
20. Add 50 μ L of Elution Buffer (65 °C Buffer PE)
21. Incubate at 65 C for five minutes

22. Centrifuge at 1 KxG for one minute
23. Repeat previous three steps with another 50 μ L of Elution Buffer and elute into the same tube
24. Freeze at -20 °C

Reagents

Red Algal Buffer for 100 mL (from Saunders 1993 JPhyc 29:251-254)

| Reagent | Amount | [Final] |
|----------------------------|----------|---------|
| 1M Trizma | 4.48 mL | 0.448M* |
| 1M Tris HCl | 5.52 mL | 0.552M* |
| 0.5 M Na ₂ EDTA | 10 mL | 0.05M |
| NaCl | 1.165 gm | 0.2M |
| Potassium Acetate | 24.54 gm | 2.5M |

* = 0.1M total Tris

Brown Algal Buffer for 100mL (from McDevit and Saunders 2009 JSP:132)

| Reagent | Amount | [Final] |
|----------------------------|----------|---------|
| 1M Tris base | 4.48 mL | 0.448M* |
| 1M Tris HCl | 5.52 mL | 0.552M* |
| 0.5 M Na ₂ EDTA | 10 mL | 0.05M |
| NaCl | 1.169 gm | 0.2M |
| CaCl ₂ | 3.329 gm | 0.3M |

* = 0.1M total Tris

Adjust to pH 8.0 and top off to 100ml

Other Buffers

A Machery-Nagel NucleoSpin Plant II kit for DNA, RNA, and protein provided the other buffers.

ProK and Tween diluted from stock